

The PDA Kiosk Software can be used to lockdown the PDA and prevent users from accessing other applications. **Administrators can tap and hold the background to unlock the PDA**, where as all other users can only access the programs listed on the screen. To install the software please visit <https://web.oqanalyst.com/PDA/>

By default the software is setup to only be able to access the OO Analyst application and is configured with the default Admin Code of **2673** that can be used to unlock the device for troubleshooting or configuring other settings such as WiFi.

The Admin Code can be changed by unlocking the device and selecting Start -> File Explorer -> Browse to “My Device” at the top -> Program Files -> PocketKiosk and selecting PKPasscode. After updating the Admin Code, insert the tip of the stylus into the hole located on the bottom or back of the PDA to restart the PDA.

Whenever the PDA restarts, a countdown timer will ask if you want to continue on to the Kiosk or unlock the device. If you select Yes or let the timer run out, the PDA will lockdown and load the Kiosk Software. If you select No you can enter the admin code to unlock the device and load into the default today screen.

If you would like to remove the Kiosk software, you can unlock the device and use the Remove Programs wizard under the System Settings to uninstall the “OO Measures PDA Kiosk” program. If you wish to add additional programs to the Kiosk home page please contact OO Measures for assistance.

